

# User Manual

## QualityMix N<sub>2</sub>O

### Oxygen-Nitrous Oxide Mixer

QualityMix N<sub>2</sub>O 70

QualityMix N<sub>2</sub>O 50

QualityMix N<sub>2</sub>O 50 FIX



*Keep these instructions!*



DEHAS Medical Systems GmbH  
Wesloer Straße 107-109  
23568 Lübeck, Germany



Tel: (+49) 451 80904-0  
Fax: (+49) 451 80904-111  
[www.DEHAS.de](http://www.DEHAS.de)








# Contents










- 1. Explanation of the most important abbreviations..... 2
- 2. Safety information - warnings, safety precautions and labelling information  
2
- 3. Package contents and inspection upon receipt..... 4
- 4. Intended application..... 4
- 5. Before first use ..... 4
- 6. Technical Data ..... 6
- 7. Pressure drop in the system..... 7
- 8. Transport and storage requirements ..... 8
- 9. Dryness and composition of gas supplies ..... 8
- 10. Illustrations and naming of components..... 8
- 11. Installation ..... 11
- 12. Alarm test..... 11
- 13. Initial operation ..... 12
- 14. Cleaning and disinfection..... 13
- 15. Maintenance..... 15
- 16. Return of Goods ..... 16
- 17. Disposal ..... 16
- 18. Troubleshooting..... 16
- 19. Warranty conditions..... 17

## 1. Explanation of the most important abbreviations

FI <sub>O</sub> <sub>2</sub>	Fractional concentration of inspiratory oxygen
DISS	Diameter Index Safety System
NIST	Non-interchangeable Screw Thread System
Bar	Unit of measurement for pressure
l/min	Litres per minute

## 2. Safety information - warnings, safety precautions and labelling information

Symbol	Description
	This symbol indicates that the device complies with the requirements of Regulation 93/42/EEC regarding medical devices and all applicable international standards.
 <b>WARNING</b>	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
 <b>ATTENTION</b>	Use of this symbol indicates a potentially hazardous situation which, if not avoided, could result in equipment damage.
 or 	Indicates the need for the user to consult the user manual.
	Follow the instructions!
	DO NOT USE OIL

<p style="text-align: center;"><b>SN</b></p>	<p>Indicates the serial number of the manufacturer so that a specific medical device can be identified.</p>
<p style="text-align: center;"><b>REF</b></p>	<p>Indicates the manufacturer's order number so that the medical device can be identified.</p>
<p style="text-align: center;"><b>MD</b></p>	<p>Medical device</p>
<p style="text-align: center;"><b>UDI</b></p>	<p>UDI Barcode</p>
<p style="text-align: center;"></p>	<p>Non-sterile</p>
<p style="text-align: center;"></p>	<p>Date of manufacture</p>
<p style="text-align: center;"></p>	<p>Indicates the manufacturer of the medical device according to EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.</p>
<p style="text-align: center;"></p>	<p>Indicates a medical device that should not be used if the packaging is damaged or opened.</p>
<p style="text-align: center;"></p>	<p>Designates a medical device that must be protected against moisture.</p>
<p style="text-align: center;"></p>	<p>Designates a medical device that may break or be damaged if handled carelessly.</p>
<p style="text-align: center;"></p>	<p>Describes the temperature limits to which the medical device can be safely exposed.</p>
<p style="text-align: center;"></p>	<p>Indicates the humidity range to which the medical device can be safely exposed.</p>
<p style="text-align: center;"></p>	<p>Indicates the range of atmospheric pressure to which the medical device can be safely exposed.</p>

### **3. Package contents and inspection upon receipt**

- Package contents:
- 1 Basic unit consisting of:
    - 1 1 mixer with adjustment unit
    - 1 Nitrous oxide module with flush function (with a permanently installed demand valve, the nitrous oxide module is not absolutely necessary, as the flush function is available in the demand valve)
  - 2 Connection hoses (O<sub>2</sub> & N<sub>2</sub>O)
  - 2 Pressure reducer 3.8 Bar (O<sub>2</sub> & N<sub>2</sub>O)
  - 1 User manual

#### **According to order option**

##### **1 Demand System N<sub>2</sub>O with connection hose**

**or**

##### **1 Flow meter with connection block and reservoir bag**

Inspection: Remove the device from its packaging and inspect it for damage. If you notice any damage, DO NOT use the device and contact your distributor.

### **4. Intended application**

The QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer is designed for the administration of a continuous and accurate mixture of medical nitrous oxide and medical oxygen to infants, children and adults via the exit port. The exact fractionated inspiratory oxygen nitrous oxide concentration (FIO<sub>2</sub>/FIN<sub>2</sub>O) corresponds to the selected FIO<sub>2</sub>/FIN<sub>2</sub>O setting on the control knob (rotary selector).

#### **Indication:**

This device is to be used by patients under the supervision of trained specialist personnel for pain therapy.

#### **Contraindication:**

Do not use with patients who cannot breathe independently. Do not use for life support or lifesaving.

### **5. Before first use**

<b><i>Read all instructions before use!</i></b>
---

These instructions for use provide qualified personnel with instructions for installing and operating the QualityMix N<sub>2</sub>O. The instructions are intended for your safety and to protect the device from damage. If you do not understand any

information or instructions in this instruction manual, do not use the device and contact your supplier.

 **DANGER**

This product is not intended for use as a life-saving or life-supporting device.

 **WARNING**

- The QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer should only be operated by medical personnel under the direct supervision of a licensed physician
- The QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer should be used only for the purposes described in these operating instructions.
- Review the prescribed dose before administration to the patient and monitor the administration frequently
- The QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer may **only** be serviced by a qualified technician
- Always comply with EN and DIN standards for the handling of medical gas products, flow meters and oxygen
- The oxygen/nitrous oxide concentration must be confirmed with an oxygen analysis/monitoring device
- **DO NOT** interfere with the alarm
- **DO NOT** use the mixer when the alarm sounds
- **DO NOT** use oil in or near the mixer
- **DO NOT** use the mixer near flames, combustible/explosive materials, vapours or gases
- **NEVER** smoke in an area where oxygen is administered
- When administering O<sub>2</sub>/N<sub>2</sub>O, personnel may be exposed to N<sub>2</sub>O. The applicable national laws specifying exposure limits and safety regulations for handling medical gases in the workplace must be complied with. Possible exposure can be prevented by continuous, effective control of the system, ventilation and working practices
- The oxygen concentration rotary switch cannot be rotated 360 degrees. Turning the switch to less than 30% or more than 100% oxygen will damage the mixer

 **WARNING**

- Close the gas supply when the QualityMix N<sub>2</sub>O oxygen/ nitrous oxide mixer is not being used.
- Always ensure an adequate supply of the gases used. After each use, check the level of the gas source before reuse.

- Store the QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer in a clean, dry place when not in use.
- The QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer does NOT contain magnetic materials containing iron and is MRI compatible (max 3 Tesla). A clearance of 2 meters must be maintained.

**The O<sub>2</sub> monitor as an accessory is not suitable for MRT.**



**ATTENTION**

- Always ensure that all connections are **secure** and **tight**
- **NOT** suitable for sterilization
- **DO NOT** immerse in liquids
- **DO NOT** sterilize with ethylene trioxide (EtO)
- **DO NOT** use if dirt or impurities are present on or near the mixer or connectors
- **DO NOT** clean with aromatic hydrocarbons
- The inlet pressure of the supply source must correspond with the specifications of the mixer.
- **When using any gas source, always use a pressure reducer within a pressure range of 3.2 to 3.8 bar for nitrous oxide (N<sub>2</sub>O) and for oxygen (O<sub>2</sub>) with a pressure of 3.6 to 4.2 bar. A pressure difference of approx. 0.4 bar higher O<sub>2</sub> pressure is required. (The pressure reducers are included in the scope of delivery for gas terminals in hospitals)**

**6. Technical Data**

<b>Model:</b>	<b>QualityMix N<sub>2</sub>O Variants:</b> <b>QualityMix N<sub>2</sub>O 50 FIX</b> <b>QualityMix N<sub>2</sub>O 50</b> <b>QualityMix N<sub>2</sub>O 70 (Mechanical stop at 50% N<sub>2</sub>O)</b>
Main output flow	1-100 l/min minimum, depending on the primary pressure
Emergency flow (failure of nitrous oxide or oxygen supply)	>50 l/min
FLUSH Button	100% O <sub>2</sub> approx. 50 Litre/Minute

Alarm activation when supply pressure drops	Alarm on - at a pressure difference between the two gases of approx. 0.9 - 1.8 bar. Alarm off - before the pressure difference between the two gases falls below 0.7 bar.
Alarm volume	≥80 dB at a distance of 1 m.
Adjustment range of N <sub>2</sub> O concentration	0-50% or 0-70% or FIX 50% N <sub>2</sub> O depending on the model With a mechanical barrier at 50% nitrous oxide concentration in the 0-70% model. Results in an N <sub>2</sub> O value equivalent to an O <sub>2</sub> concentration of between 100% - 30% O <sub>2</sub>
Gas inlet pressure with required pressure reducer	3.2 - 3.8 bar nitrous oxide and 3.6 - 4.2 bar oxygen within max. 0.4 bar pressure differential
Accuracy of the mixed gas (FIO <sub>2</sub> )*	± 3 % oxygen
Connection types	1 x DISS output for mixed gas and 1 x NIST input for nitrous oxide N <sub>2</sub> O and oxygen O <sub>2</sub> respectively
Dimensions LxBxH	13 x 16.5 x 18.2 cm
Weight	2100 g
Operating temperature	+5°C to +50°C

The QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer has been degreased for oxygen utilization prior to delivery. The reversed gas flow of the oxygen nitrous oxide mixer complies with Clause 9 of ISO 11195:2018. The oxygen analyser used must comply with ISO 80601-2-55 and CE regulations.

## 7. Pressure drop in the system

Low flow	≤0.14 bar with inlet pressures of 3.8 bar and a flow rate of 10 l/min at >50% F <sub>IN<sub>2</sub>O</sub>
High flow	≤0.21 bar at inlet pressures of 3.8 bar and a flow rate of 30 l/min at 50% F <sub>IN<sub>2</sub>O</sub>



## 8. Transport and storage requirements

Temperature range	-20 °C to 50 °C
Humidity	max. 95% non-condensing humidity

## 9. Dryness and composition of gas supplies


### Nitrous oxide (N<sub>2</sub>O):


The medical nitrous oxide must meet all of the requirements for medical nitrous oxide (N<sub>2</sub>O) according to the European Pharmacopoeia.

### Oxygen (O<sub>2</sub>):

The oxygen used must meet all of the requirements for medical oxygen (O<sub>2</sub>) according to the European Pharmacopoeia.


## 10. Illustrations and naming of components

 <b>ATTENTION</b>
If you follow the instructions regarding processing, the labelling on the devices will remain intact. If the labelling should nevertheless become illegible or missing, please contact the manufacturer or your local contact person.





*This illustration depicts the QualityMix N<sub>2</sub>O with connected flow meter*

Component	Description
Rotary switch for oxygen and nitrous oxide concentrations	A rotary switch for setting the oxygen concentration to between 0 and 50 (70)%. The FIO <sub>2</sub> scale is for reference purposes only. This rotary switch cannot be rotated 360°. The rotary switch starts at 0% and extends to 50 (70)%.
Main outlet	A male threaded DISS connector with shut-off valve that provides gas flow when connected to a control device such as a flow meter with reservoir bag or an N <sub>2</sub> O demand system.
Connection for medical oxygen	A NIST oxygen connection with female thread and one-way valve for connecting an oxygen supply hose.
Connection for medical nitrous oxide	A male threaded NIST connector and a one-way valve to connect a nitrous oxide supply hose.
Alarm	An audible alarm that sounds in the event of excessive pressure drop or failure of oxygen and/or nitrous oxide supply.
Flush O <sub>2</sub> button	For administration of 100% O <sub>2</sub> , independent of the concentration setting on the rotary switch of the QualityMix N <sub>2</sub> O mixer.



*This figure shows the optional N<sub>2</sub>O demand system with scavenging system for the exhaled air - different versions are possible here.*

## 11. Installation

 <b>Warning</b>
<ul style="list-style-type: none"><li>• <b>Read</b> the user manual before installing or using the unit.</li><li>• <b>Monitor</b> the oxygen and nitrous oxide concentration with an oxygen analyzer/monitoring device.</li></ul>
 <b>Attention</b>
Check the QualityMix N <sub>2</sub> O mixer for visible damage before use and do <b>not</b> use if it is damaged.

**Note: Perform the following test before operating the device for the first time**

### **Alarm test (see section 12)**

The alarm test of the QualityMix N<sub>2</sub>O blender must be performed once a month.

#### **Preparation for the alarm test**

1. Attach the QualityMix N<sub>2</sub>O oxygen/nitrous oxide mixer to a rail or support rod in an upright position.
2. Connect the nitrous oxide and oxygen supply lines to the appropriate inlet ports at the bottom of the mixer.
3. Connect a flow meter or N<sub>2</sub>O demand valve to one of the output ports.

#### **Flow capacity of the output:**

- N<sub>2</sub>O mixer, all variants: 0 – 100 l/min minimum, depending on the primary pressure
4. Connect an exhaust line to the outlet port of the flow meter or to the intended connection of the N<sub>2</sub>O demand valve.

## 12. Alarm test

1. Connect the QualityMix N<sub>2</sub>O mixer to the nitrous oxide and oxygen sources. Connect the nitrous oxide source first.
2. Pressurize the mixer and turn the flow meter anti-clockwise.

3. Set the rotary switch for nitrous oxide concentration to 50% ( $\text{N}_2\text{O}$ ).
4. Hold down the flush button on the demand valve and / or open the flow meter to 5 -10 liters flow.
5. Disconnect or switch off the nitrous oxide supply to the QualityMix  $\text{N}_2\text{O}$  oxygen/nitrous oxide mixer and press the flush button on the demand valve. The mixer should make a loud beeping sound as an alarm. This sound indicates that the alarm is working properly.
6. Reconnect and activate the nitrous oxide supply to the mixer; the beeping should stop.
7. Disconnect or switch off the oxygen supply to the mixer and press the flush button on the demand valve. The beeping sound indicates that the alarm is functioning properly.
8. Reconnect and activate the oxygen supply to the mixer; the beeping sound should stop.
9. Finally turn the dial to 100%  $\text{O}_2$  and press the flush button.
10. If the alarm does not function correctly, **DO NOT USE** the unit.

### 13. Initial operation




#### Attention


Before use, check the QualityMix  $\text{N}_2\text{O}$  oxygen/nitrous oxide mixer for visible damage and do not use if it is damaged.

1. Attach the mixer to the rail or stand bracket.
2. Connect the nitrous oxide and oxygen supply lines to the mixer and supply.
3. Connect the flow meter or the  $\text{N}_2\text{O}$  demand system to the output of the mixer.
4. Set the nitrous oxide concentration rotary switch to the prescribed value.
5. Activate the  $\text{N}_2\text{O}$  gas supply (alarm sounds)
6. Activate the  $\text{O}_2$  gas supply (alarm off)

7. Check the flow of the oxygen/nitrous oxide mixture to the patient by pressing the flush button on the demand valve.
8. Check the oxygen/nitrous oxide concentration with an oxygen analysis/monitoring device.
9. When the oxygen/nitrous oxide mixer is not in use, shut off the gas supply or disconnect the appliance from the gas supply.

## 14. Cleaning and disinfection

 <b>ATTENTION</b>
<ul style="list-style-type: none"> <li>• <b>NOT suitable for sterilization.</b></li> <li>• <b>NEVER immerse the QualityMix N2O oxygen-nitrous oxide mixer in liquids.</b></li> <li>• <b>DO NOT use strong solvents or abrasive cleaners.</b></li> <li>• <b>DO NOT clean with aromatic hydrocarbons.</b></li> </ul>

<p>Manufacturer:</p> <p>DEHAS Medical Systems GmbH</p> <p>Wesloer Straße 107-109</p> <p>23568 Lübeck</p>	<p>Sterilization method:</p> <p>N/A</p> <p>The product is not intended for the sterilization process.</p>	
--	---	---

### Product described:

QualityMix N2O Serie

<b>WARNINGS:</b>	<p>Do not use disinfectants containing phenol.</p> <p>Do not use strong solvents or abrasive cleaners.</p> <p>Do not clean with aromatic hydrocarbons.</p> <p>Do not autoclave!</p> <p>Do not sterilize!</p> <p>Do not immerse in liquids!</p>
------------------	--

The instructions listed have been validated as SUITABLE by the medical device manufacturer for preparing a medical device for reuse. Standard: ISO 17664-2 Reprocessing of health care products - Information to be provided by the medical device manufacturer for the reprocessing of medical devices - Part 2: Non-critical medical devices.

The reprocessor is responsible for ensuring that the reprocessing actually carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired results.

<b>Instructions</b>	
	The outside of the unit must be wiped clean at regular intervals or at the latest after each patient in accordance with the applicable hygiene standards.
<b>Preparation for cleaning and disinfection:</b>	All gas connections from the device must be disconnected before cleaning and disinfection.
<b>Cleaning: Manuel</b>	For this procedure: <ol style="list-style-type: none"> <li>1. Wipe the outer surfaces with a lint-free cloth moistened with a mild detergent and water.</li> <li>2. Wipe dry with a dry, lint-free cloth.</li> </ol>
<b>Disinfection: Manuel</b>	<ol style="list-style-type: none"> <li>3. Wipe the surface of the product with a lint-free cloth moistened with non-oxidizing disinfectant.</li> <li>4. Observe the exposure time of the disinfectant manufacturer according to the required spectrum of activity.</li> </ol>
<b>Drying:</b>	<ol style="list-style-type: none"> <li>5. After the exposure time specified by the disinfectant manufacturer, wipe dry with a dry, lint-free, low-germ cloth.</li> <li>6. Ensure that the product is free of disinfectant residues.</li> </ol>
<b>Manufacturer's Recommendation:</b>	The manufacturer recommends the use of the disinfectant Bacillol <sup>®</sup> 30 Foam, Bacillol <sup>®</sup> 30 Tissues, Bode Chemie GmbH & Co. The current product data sheet of the disinfectant manufacturer must be observed.php
<b>Maintenance, control and testing:</b>	After each cleaning and disinfection, a visual and functional check must be performed (see chapter 13).
<b>Packaging:</b>	The device must be kept in the packaging validated by DEHAS for transport.
<b>Storage:</b>	<p>Store the product in a dry, clean place free from lubricants, oil and other contaminants.</p> <p>Storage temperature should not fall below or exceed -20°C to 50°C.</p> <p>Humidity max. 95% non-condensing.</p>
<b>Transport:</b>	N/A
<b>Additional information:</b>	N/A

<b>Contact:</b>	DEHAS Medical Systems GmbH, Wesloer Straße 107-109, 23568 Lübeck, Germany  Tel.: +49 451 80904-0, Fax: +49 451 80904-111, E-Mail: info@dehas.de
-----------------	---

## 15. Maintenance

**The following maintenance and inspection tasks must be carried out:**

- **Monthly** check of the alarm by the user
- A **safety inspection** (SI) must be carried out **every** year by a **trained operator** or medical technician.
- Maintenance should be carried out **at least every 2 years** by trained specialist personnel. Inspection of the reversed gas flow is part of the maintenance and is therefore carried out every 2 years.

### **Reverse gas flow test**

1. Set the nitrous oxide concentration of the mixer to 50%.
2. Connect the nitrous oxide supply hose to the mixer and the gas supply and open the supply.

Measure the flow rate at the oxygen inlet using a suitable measuring instrument.

The flow rate should not exceed 10 ml/h.

If the flow rate is greater than 10 ml/h, the duckbill valve in the oxygen inlet must be replaced in accordance with the service instructions and the measurement repeated.

3. Connect the oxygen supply hose to the mixer and the gas supply, and open the supply.

Measure the flow rate at the nitrous oxide inlet using a suitable measuring instrument.

The flow rate should not exceed 10 ml/h.

If the flow rate is greater than 10 ml/h, the duckbill valve in the nitrous oxide inlet must be replaced in accordance with the service instructions and the measurement repeated.



## 16. Return of Goods

Please contact your distributor regarding this matter. They will coordinate the return for you. It is important that you provide a description of the problem so that the return can be processed in a targeted manner. All returns must be sent in sealed containers to avoid damage. The distributor is not responsible for equipment that is damaged during transport.

## 17. Disposal

This device and its packaging do not contain any hazardous substances. No special precautions are required when disposing of the device and/or its packaging.

**Please recycle**

## 18. Troubleshooting

If the oxygen/nitrous oxide mixer fails, refer to the troubleshooting section below. If this does not solve the problem, please contact your local distributor.

<b>Problem</b>	<b>Possible cause</b>	<b>Remedy</b>
<b>Discrepancy between the setting of the oxygen concentration on the mixer and on the analysis/monitoring device (more than 3 %)</b>	1. Pre-pressure too unequal/too low	Pre-pressure inspection: optimum pre-pressure 3.9 to 6.5 Bar
	2. Analysis/monitoring device does not accurately register	Recalibrate the monitoring device or use a different analysis/calibration monitoring device
	3. Gas supply contaminated	Check gas supply with calibrated oxygen analyser/monitoring device to ensure oxygen content is 100% and nitrous oxide content is 0%
	4. The Flow of downward mounted device causes backflow or restricted flow	Disconnect the mixer. Check the oxygen concentration at the mixer outlet

<b>No flow at the mixer output</b>	1. Gas supply turned off	Turn on gas supply
	2. Gas supply not connected	Connect gas supply
<b>Alarm sounds</b>	Difference between oxygen and nitrous oxide inlet pressures higher than prescribed	Correct pressure difference until oxygen and nitrous oxide pressures meet specifications

**19. Warranty conditions**

The distributor guarantees that the mixer is free from defects in workmanship and/or materials for the following period of time:

One (1) year after delivery

The distributor will, with written notice and with evidence that the device has been stored, installed, serviced and operated in accordance with instructions and standard industry practices and that no modifications, substitutions or modifications have been made to the product, will correct such defect by appropriate repair or replacement at the distributor’s expense.

VERBAL STATEMENTS DO NOT CONSTITUTE A GUARANTEE.

The distributor is not authorized to provide oral guarantees about the product described in this agreement, and such statements are not binding and are not part of the purchase agreement. Therefore, this second statement is the final, complete and exclusive representation of the terms of this agreement.

The current version of the General Terms and Conditions of the distributor and German law shall apply.

# Declaration of Conformity



DEHAS Medical Systems GmbH  
Wesloer Straße 107-109  
23568 Lübeck  
GERMANY



0482



QualityMix N<sub>2</sub>O 50 FIX

QualityMix N<sub>2</sub>O 50

QualityMix N<sub>2</sub>O 70

And the relevant accessories

**Classification:** IIb

**Classification Criteria:** Clause 3.2 Rule 11 of Annex IX of the MDD

We hereby declare under our sole responsibility that the products mentioned above comply with the provisions of the following directives and standards of the EC Council. All documents are retained on the premises of the manufacturer and the notified body.

**Guidelines:** General application guidelines: Medical Device Directive (MDD), Council Directive 93/42/EEC of 14 June 1993 Annex II.3 on medical devices of the European Parliament.

<b>Applied standards:</b>	DIN EN 1041	ISO 11195
	EN ISO 14971	ISO 18562-1
	DIN EN ISO 15001	ISO 18562-2
	DIN EN ISO 15002	ISO 18562-3
	DIN EN ISO 15223-1	ISO 10993-1
	DIN EN 62366-1	

**Notified Body:** DNV Medcert GmbH /  0482

**Address:** Pilatuspool 2, 20355 Hamburg; GERMANY

**Certificate number:** 4153DE410200327      Expiry date: 05/2024

**Already manufactured devices:** Traceability via serial number

**Valid from/to:** 27-03-2020 to expiry date

**Manufacturer representative:** Jens Mittendorf

**Position:** CEO /Head of Development

**Date of issue:** 03-04-2020

**Your contact person for sales and service:**

A large, empty rectangular box with a thin black border, intended for the user to provide contact information for sales and service.