

Budget Mix

Low-Flow O₂-Air-Blender
 High-Flow O₂-Air-Blender (fig.)



BudgetMix - The compact and affordable O₂-AIR Blender from DEHAS. Reduced to the essentials, but without loss of precision. A 2-chamber mixing system for oxygen and compressed air that does what you expect. Like its big brother - the BudgetMix regulates the set value between 21% and 100% oxygen content precisely and quickly - but without any further upgrade options

Available accessories: Flowmeter 0-3/0-6/0-15/3.2-32/8.5-85 lpm, silencer, rail mount.

Value	Budget Mix Low-Flow	Budget Mix High-Flow
Measurements (WxHxD – without connectors)	63 x 115 x 115 mm	63 x 115 x 115 mm
Weight Base Unit	ca. 1400 g	ca. 1400 g
Mixing range	21 - 100 % O ₂	21 - 100 % O ₂
Output Flow	0-27 l/min	0-105 l/min
Connection Inlet	NIST	NIST
Connection Outlet	G 9/16" DISS	G 9/16" DISS
Accuracy	+/- 3 % Sauerstoff	+/- 3 % Sauerstoff
Volume of the alarm	≥80 dB at a distance of 30 cm	

CE-Sign/Cert.number: 13485:2016 certified

CE 0482

WE THINK IN PRESSURE & FLOW