

Flow Meter directly pluggable

for flows from 3 LPM up to 85 LPM



We supply our directly pluggable flow meters in a wide range, both in terms of the available flow rates and the possible international standards.

All flow meters are available for oxygen and air.

Security features:

Device type: not electrically operated
Function: mechanical / pneumatic
CE mark / number: 13485: 2016 certified



Available standards:

German (DIN), French (NF), British (BS), Scandinavian (SS, AGA), Switzerland (CH, Carbamed), Italy (UNI).

| Available flow rates | Available gas types | Nominal pressure | Deviation | Outlet |
|----------------------|----------------------|------------------|-----------|------------|
| 0 - 3 LPM | O ₂ , AIR | 4,5 bar | < 10% | DISS 9/16" |
| 0 - 6 LPM | O ₂ , AIR | 4,5 bar | < 10% | DISS 9/16" |
| 0 - 15 LPM | O ₂ , AIR | 4,5 bar | < 10% | DISS 9/16" |
| 3 - 32 LPM | O ₂ , AIR | 4,5 bar | < 10% | DISS 9/16" |
| 8 - 85 LPM | O ₂ , AIR | 4,5 bar | < 10% | DISS 9/16" |

WE THINK IN PRESSURE & FLOW