


Gas supply hoses

Standard goods and custom-made products



We have been manufacturing connection hoses for medical patient care in our company since 2002. All hoses are manufactured as soon as the order is received and are then available at short notice. Our extensive standard program is supplemented by customer-specific special constructions. Due to our extensive parts warehouse, we can also provide special constructions at short-term.



Safety features for patient and user:

Device type: not electrically operated
Function: mechanical
Industry standard: DIN EN ISO 11195: 2018
CE mark / number: 13485: 2016 certified



Hose material:

Medical gas hose according to EN ISO 5359. Five-layer hose made of soft PVC with fabric insert
Made of high strength polyester. Antistatic and heat resistant.

Inner diameter	Outer diameter	Thickness	Pressure resistance	Weight
6,3 mm	12,4 mm	3,0 mm	15 bar	106 g/m

Available standards:

German (DIN), French (NF), British (BS), Scandinavian (SS, AGA), Switzerland (CH, Carbamed), Italy (UNI).

Available gas types:

Oxygen (O₂), air (AIR), nitrous oxide (N₂O), carbon dioxide (CO₂), vacuum (VAC), compressed air (AIR8)

WE THINK IN PRESSURE & FLOW