


Flow meter for vertical mounting

for flows from 3 LPM up to 85 LPM



Our oxygen flow meter for mounting on mobile supply units such as infusion stands. It is suitable for pipe diameters up to 40 mm.

Mounting on the standard rail is also possible by turning the clamp by 90 °.

Security features:

Device type: not electrically operated
 Function: mechanical / pneumatic
 CE mark / number: 13485: 2016 certified



Adjustable flow	Part-no.	Nominal pressure	Deviation	Inlet	Outlet
0 - 3 LPM	D-FL-O2-3-ST-NIST	4,5 bar	< 10 %	NIST	DISS 9/16"
0 - 6 LPM	D-FL-O2-6-ST-NIST	4,5 bar	< 10 %	NIST	DISS 9/16"
0 - 15 LPM	D-FL-O2-15-ST-NIST	4,5 bar	< 10 %	NIST	DISS 9/16"
3 - 32 LPM	D-FL-O2-32-ST-NIST	4,5 bar	< 10 %	NIST	DISS 9/16"
8 - 85 LPM	D-FL-O2-85-ST-NIST	4,5 bar	< 10 %	NIST	DISS 9/16"

WE THINK IN PRESSURE & FLOW