

QualityREG Cylinder pressure regulators

Regulators for medical gases



In addition to the standard program with stepped or stepless pressure reducers, we offer pressure regulators with customer-specific connections according to your requirements. Our focus is on the supply of oxygen (O₂), but other medical gases such as AIR / CO₂ / N₂O are also possible.

We offer pressure reducers for medical gases that can be configured according to your needs. Contact us!



Technical Data	Value
Max. pressure on the input side	200 bar
Max. pressure on the output side	See type label
Max. outlet pressure of the separate pressure outlets:	4 bar – 5 bar
Flow accuracy	Either ±10 % or ±0,5 l/min (the higher value is valid)
Screw thread cylinder adapter	Similar to the type of gas and national standard
Connection constant outlet	See list of variations
Screw thread variable outlet	9/16" UNF
Caliber hose connector adapter	6 mm according to DIN-EN-ISO 13544
Weight ; depending on type	560 g to 830 g
Allowed temperature range	–20 – +60 °C
Type of gas	See type label

CE 0482

WE THINK IN PRESSURE & FLOW