

QualityREG INLINE, adjustable

Adjustable tapping point pressure reducer



Our solution for the operation of compressed air driven devices or for the compensation of very different line pressures for the operation of gas mixers or the compensation of fluctuating line pressures.

The pressure reducer is available for all non-corrosive gases. The device-side connection is made via a NIST screw connection.

Security features:

Device type: not electrically operated
Function: mechanical / pneumatic
Industry standard: DIN 13260 part 2
CE mark / number: 13485: 2016 certified



Technical specifications	Value
Inlet pressure	max. 8 bar
Outlet pressure	continuously adjustable 0,1 - 6 bar
Manometer	0 - 6 bar
Deliverable standards	German (DIN), French (NF), British (BS), Scandinavian (SS, AGA), Switzerland (CH, Carbamed), Italy (UNI).

WE THINK IN PRESSURE & FLOW